




Closed Line Adaptor

REF CLA100



Instructions for Use

Description

The Halo® Closed Line Adaptor is a needle-free, luer lock port adaptor for the closed transfer of antineoplastic and hazardous drugs from vial to patient in a healthcare setting. Use required and accepted IV practice. Do not reuse to avoid contamination. Sterile (ethylene oxide). Sterile unless package is damaged. Rx Only.


Indication for Use

The Halo® system is an airtight and leak proof closed system drug transfer device (CSTD) that mechanically prohibits the transfer of environmental contaminants into the system and the escape of drug or vapor concentrations outside the system, thereby minimizing individual and environmental exposure to drug vapor, aerosols and spills. The Halo® system also prevents microbial ingress for up to 7 days.

Precaution

Always grip white components when using the Halo® System. Do not use Halo® Closed Line Adaptor as protection cap during transport or storage.

The sealing membranes are validated for up to 14 perforations. The performance of the sealing membrane may be reduced after multiple perforations. Discard after 14 perforations.



- Instructions**
1. Attach the Closed Line Adaptor to the patient's luer-locking IV port or point of access.
 2. Hold and Push the Closed Syringe Adaptor onto the connection interface of the Closed Line Adaptor until fully engaged.
 3. Perform drug transfer
 4. Hold and Pull the Closed Syringe Adaptor to separate from the connection interface.
 5. Dispose without disconnecting the Halo® components.

Symbol Legend



Catalogue number



Batch code



Do not re-use



Do not use if package is damaged



Temperature limits to which the device can be safely exposed



Caution, consult accompanying documents



Sterilized using ethylene oxide



Use by date



Caution: Federal law (US) restricts this device to sale by or on the order of a physician



Manufacturer
Corvida Medical, Inc.



Consult instructions for use



Non-pyrogenic fluid path



Not made with natural rubber latex



Does not contain the Phthalate
DEHP

Distributed By



Fresenius Kabi USA, LLC
Three Corporate Drive
Lake Zurich, Illinois 60047
www.fresenius-kabi.us
+1 800-933-6925